

# IN BED PATIENT SAFETY PRODUCTS

See reverse for Entrapment Zones Diagram

## WEDGE SIDE RAIL PADS (HEAVY DUTY)

Entrapment Zone: 1, 2, 3 & 5  
 Item #8600-HDF - Full Length: 14"H x 58"L - c/w Velcro Straps  
 Item #8600-HDH - Half Length: 14"H x 27"L - c/w Velcro Straps



- Designed specifically to reduce patient entrapment and unintended bed exits in several high risk (1, 2, 3 & 5) entrapment zones when using side rails
- Provides added padding to side rails to reduce patient injuries
- Can be used for patient side positioning
- Eliminates the need for raised perimeter mattresses
- Wipable Recovery 5 cover

## WEDGE SIDE RAIL PADS (STANDARD)

Entrapment Zone: 1, 2 & 3  
 Item #8600F - Full Length: 17"H x 58"L - c/w Velcro Straps  
 Item #8600H - Half Length: 17"H x 27"L - c/w Velcro Straps



- Designed specifically to reduce patient entrapment and unintended bed exits in several high risk (1, 2 & 3) entrapment zones when using side rails
- Eliminates the need for raised perimeter mattresses
- Wipable Recovery 5 cover

## CUSTOM BED RAIL GAP FILLER

Entrapment Zone: 1  
 Item #8602  
 Size: Call for a custom size quotation



- Designed specifically to prevent patient entrapment in Zone 1
- Prevents patient movement through side rails
- Wipable Recovery 5 cover
- Custom made to fit most side rails

## IN BED PATIENT POSITIONING AND SAFETY BOLSTERS


Entrapment Zone: 1, 2, 3, 4 & 5  
 Item #8603F - Full Length: 2" x 14"H x 58"L  
 Item #8603H - Half Length: 2" x 14"H x 27"L



- Designed specifically for when bed rails are not permitted or used
- Used for the prevention of patient entrapment and unintended bed exits in high risk (1, 2, 3, 4 & 5) entrapment zones
- Eliminates the need for raised perimeter mattresses
- Non-slip base along with the patient's weight keeps the system securely in place
- Side bolsters can be used for patient side positioning and to encourage the patient to remain centered
- System secured with heavy duty straps and sturdy buckles
- Wipable Recovery 5 Cover
- 20" of sleeping surface space between each safety bolster

## VITA 5 RAISED PERIMETER MATTRESS

Entrapment Zone: 1, 2 & 3  
 Item #3050



- The Vita 5 was designed to encourage patients to position themselves in the mattress center and discourage rolling or unattended bed exits
- The raised perimeter edge reduces the gaps between the bedrail and the mattress, limiting the risk of patient impingement of extremities
- Centre cut-outs in the sturdy side bolsters make it easy to facilitate transfer and patient care
- Wipable Recovery 5 Cover

## ZONE 6 AND 7 PATIENT PROTECTION SYSTEM

Entrapment Zone: 6 & 7  
 Item #8601 - 10.5"H x 36"W x 14"L



- Helps reduce the risk of patient entrapment between the mattress and headboard and/or footboard
- High density foam protects patients from injury
- System secured with heavy duty straps and sturdy buckle
- Wipable Recovery 5 cover

## FULL LENGTH LATERAL POSITIONER

Entrapment Zones: 1 & 5  
 Item #6279 Wedge - 8" x 8" x 58"  
 Item #6276 Half Moon - 4" x 8" x 58"



- Designed to provide comfortable patient positioning while offering the patient a gentle reminder they are nearing the outside mattress edge
- Provides protective padding to side rails to reduce patient injury
- Replaces the need for bed rail padding
- Non-Slip base helps secure bolster into place while positioning the patient
- Optional straps to secure to side rails
- Wipable Recovery 5 cover

## BED RAIL PROTECTION PADS

Entrapment Zone: 1  
 Item #6277 - 4" x 8" x 27" - c/w 3 Velcro Straps  
 Item #6277A - 4" x 8" x 34" - c/w 3 Velcro Straps  
 Item #6277L - 4" x 8" x 72" - c/w 4 Velcro Straps



- Half moon design provides the patient with a gentle reminder they are approaching the outside mattress edge
- High quality firm foam provides protection for patients from hard bed rails
- Adjustable straps with Velcro offers flexible fitting capabilities for most bed rail configurations
- Non-Slip material reduces pad movement
- Wipable Recovery 5 cover

## MATTRESS EXTENSION BOLSTERS

Entrapment Zones: 2, 3, 4 & 7  
 Item #6278 (Custom Bolster)  
 Size: Call for a custom size quotation



- Designed to correct improperly fitted mattresses and hospitals beds to reduce the possibility of patient entrapment
- Eliminates the gap between the mattress and the headboard, footboard or side rails
- Eliminates the need to use towels, blankets or pillows
- Non-Slip base keeps the bolster firmly in place
- Wipable Recovery 5 cover

[www.vitacaremedical.com](http://www.vitacaremedical.com)

Toll Free Phone: 1-800-263-9068 / Toll Free Fax: 1-866-622-0986

\*Most of our products can be modified to meet unique sizing requirements.