


Important: There is no weight limit when combining a minimum of 3 and 4 Star mattress inserts with Vitacare V4's, V4 RE-ACT's or leveling pads. 1 and 2 Star mattress inserts combined with any V4, V4 RE-ACT or leveling pads offers premium pressure distribution with a weight limit of 250LBS. We highly recommend for patients that spend more than 12-14 hours per day in bed use a minimum of 2 Star mattress inserts to not only provide better skin protection but also to prolong the integrity of the foam components. Rotating the upper foam sections is recommended to extend the life and performance of the foam sections.

Vitacare V4 Mattress Systems

35" x 80"



V4 Insert Options—Recovery 5 Covers

ITEM	SIZE
Vita 1935	Foam Insert 19" x 35" x 3.25"
Star 1934	Star Mattress Insert 19" x 34" x 4"
Recovery 5 Cover Options	
VC-COVER9	Elastic Bottom Mattress Cover (Shower Cap)
VS-COVER1934	Star Mattress Insert Cover 19" x 34" x 4" with 270° Zipper
VS-COVER1934 DO	Star Mattress Insert Cover 19" x 34" x 4" with double opening

Note:
To add extra longevity and ease of cleaning of the mattress and the Star mattress insert we recommend adding the protection of a Recovery 5, MRSA Platinum, or a Heavy Duty LTC cover.
Vitacare can *customize* any V4, V4 RE-ACT Mattress System i.e. double, twin, queen, king to meet your patient's and their companions needs i.e. Partners in the same bed.

V4 Insert Options—Platinum MR Covers

ITEM	SIZE
Vita 1935 MR	Foam Insert 19" x 35" x 3.25"
Star 1934	Star Mattress Insert 19" x 34" x 4"
Platinum MR Cover Options	
VC-COVER9-MR	Elastic Bottom Mattress Cover (Shower Cap)
VS-COVER1934-MR	Star Mattress Insert Platinum MR Cover 19" x 34" x 4" with 270° Zipper
VS-COVER1934DO-MR	Star Mattress Insert Platinum MR Cover 19" x 34" x 4" with double opening


Cover Option: MRSA Resistant Platinum Fabric

Specially treated Platinum MR fabric utilizes a powerful antimicrobial agent able to fight the most aggressive bacteria and fungi (MRSA) to restrict the growth of the MRSA strains. Stretchy incontinent fabric reduces skin shear and friction. Odor is controlled by preventing bacterial and fungal growth.

NOTE: To order your V4 Products and accessories with MR Platinum Fabric use the MR product Item numbers listed below and lower right side.



Star Mattress Insert 19" x 34" x 4"
The Star Mattress Inserts come standard in a 4" Cell height. This allows added immersion, for increased pressure redistribution and enhanced skin protection.



Periodically rotate all top foam sections for enhanced performance and longevity

*Upgraded Platinum MRSA Cover

NOTE: The mattress cover will be adversely affected by leaving the cleaning disinfectants on the fabric for a prolonged duration or a failure to adequately rinse the agent from the fabric after recommended dwell time.

