



NOTE:
 Star 1934= Star Mattress Insert
 VLP = Vitacare Leveling Pad

Vitacare Leveling Pad Overlay Configurations

34" x 80"




1 x Star 1934
 3 x VLP2134




2 x Star 1934
 2 x VLP2134

ITEM	SIZE
VLP2134	21" x 34"
VLP3134	31" x 34"
VLP3480	34" x 80"
Star 1934	19" x 34" x 4"
Recovery 5 Cover Options	
VC-COVER4	Elastic Bottom Cover 36" x 80" x 4"
VC-COVER9	Elastic Bottom Cover 36" x 80" x 9"
VS-COVER 1934	Star Mattress Insert Cover 19" x 34" x 4" with 270° Zipper
VS-COVER 1934 DO	Star Mattress Insert Cover 19" x 34" x 4" with Double Opening



3 x Star 1934
 1 x VLP2134



1 x Star 1934
 1 x VLP3480



4 x Star 1934



1 x Star 1934
 2 x VLP3134

Important: There is no weight limit when combining 3 or 4 Star air cell mattress inserts combined with leveling pads to make an overlay. 1 to 2 Star inserts combined with leveling pads offers premium pressure distribution with a weight limit of 250lbs. We suggest patients that spend more than 12-14 hours per day in bed use 2 or more Star Mattress Inserts to not only provide better Skin protection but also to prolong the life of the foam components.

