




# Huntington's Patient In-Bed Protection System


- Provides patient with a safe environment
- Anti-Microbial, easy to clean nylon cover
- Drop down side bolsters for ease of patient transfers / exits
- Custom built to fit any size bed
- System sits on top of existing mattress and is secured to side rails, bed frame and the head / foot boards
- Easy to install

Premium Patient Protection with High Density Foam



Collapsible Side Bolster for ease of Patient Exiting and ADL's

Optional Wedge Shape Bolster Design



Side Rail Straps secure to Bed Rails

**Note:** MRSA resistant fabric is available